

# PRODUCT REFERENCE GUIDE

DRUGS OF ABUSE									
Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
AMP Amphetamine	Urine	Strip	DAM-101	✓	✓		1000 ng/mL	50	4.0 mm
			DAM-U101	✓			500 ng/mL	50	
			DAM-A101	✓	✓		300 ng/mL	50	
			DAM-111	✓			1000 ng/mL	25 Tests/Canister	
		Device	DAM-102	✓	✓		1000 ng/mL	40	
			DAM-A102	✓	✓		300 ng/mL		
		Dip Card	DAM-114	✓	✓		1000 ng/mL	25	4.0 mm
			DAM-U114	✓			500 ng/mL	25	
			DAM-A114	✓	✓		300 ng/mL	25	
BAR Barbiturates	Urine	Strip	DBA-101	✓	✓		300 ng/mL	50	4.0 mm
			DBA-111	✓				25 Tests/Canister	
		Device	DBA-102	✓	✓			40	
		Dip Card	DBA-114	✓	✓			25	
BUP Buprenorphine	Urine	Strip	DBU-101	✓	✓	✓	10 ng/mL	50	4.0 mm
			DBU-111	✓				25 Tests/Canister	
		Device	DBU-102	✓	✓	✓		40	
		Dip Card	DBU-114	✓	✓	✓		25	
BZO Benzodiazepines	Urine	Strip	DBZ-101	✓	✓		300 ng/mL	50	4.0 mm
			DBZ-A101	✓			200 ng/mL	50	
			DBZ-111	✓			300 ng/mL	25 Tests/Canister	
		Device	DBZ-102	✓	✓		300 ng/mL	40	
			DBZ-A102	✓			200 ng/mL		
		Dip Card	DBZ-114	✓	✓		300 ng/mL	25	
DBZ-A114	✓			200 ng/mL	25				

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
COC Cocaine	Urine	Strip	DCO-101	✓	✓		300 ng/mL	50	4.0 mm
			DCO-U101	✓	✓		150 ng/mL	50	
			DCO-111	✓			300 ng/mL	25 Tests/Canister	
		Device	DCO-102	✓	✓		300 ng/mL	40	
			DCO-U102	✓	✓		150 ng/mL		
		Dip Card	DCO-114	✓	✓		300 ng/mL	25	
DCO-U114	✓		✓		150 ng/mL	25			
KET Ketamine	Urine	Device	DKE-102	✓			1000 ng/mL	40	4.0 mm
MDMA Ecstasy	Urine	Strip	DMD-101	✓	✓	✓	500 ng/mL	50	4.0 mm
			DMD-111	✓			500 ng/mL	25 Tests/Canister	
		Device	DMD-102	✓	✓	✓	500 ng/mL	40	
		Dip Card	DMD-114			✓	500 ng/mL	25	
MET Methamphetamine	Urine	Strip	DME-101	✓	✓		1000 ng/mL	50	4.0 mm
			DME-U101	✓	✓		500 ng/mL	50	
			DME-A101	✓			300 ng/mL	50	
			DME-111	✓			1000 ng/mL	25 Tests/Canister	
		Device	DME-102	✓	✓		1000 ng/mL	40	
			DME-U102	✓	✓		500 ng/mL		
			DME-A102	✓			300 ng/mL		
		Dip Card	DME-114	✓	✓		1000 ng/mL	25	
			DME-U114	✓	✓		500 ng/mL	25	
DME-A114	✓				300 ng/mL	25			
MOP Morphine 300	Urine	Strip	DMO-101	✓	✓		300 ng/mL	50	4.0 mm
			DMO-111	✓			300 ng/mL	25 Tests/Canister	
		Device	DMO-102	✓	✓		300 ng/mL	40	
		Dip Card	DMO-114	✓	✓		300 ng/mL	25	

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
MTD Methadone	Urine	Strip	DMT-101	✓	✓		300 ng/mL	50	4.0 mm
			DMT-111	✓			300 ng/mL	25 Tests/Canister	
		Device	DMT-102	✓	✓		300 ng/mL	40	
		Dip Card	DMT-114	✓	✓		300 ng/mL	25	
OPI Opiate 2000	Urine	Strip	DOP-101	✓	✓		2000 ng/mL	50	4.0 mm
			DOP-111	✓			2000 ng/mL	25 Tests/Canister	
		Device	DOP-102	✓	✓		2000 ng/mL	40	
		Dip Card	DOP-114	✓	✓		2000 ng/mL	25	
OXY Oxycodone	Urine	Strip	DOX-101	✓	✓	✓	100 ng/mL	50	4.0 mm
			DOX-111	✓			100 ng/mL	25 Tests/Canister	
		Device	DOX-102	✓	✓	✓	100 ng/mL	40	
		Dip Card	DOX-114			✓	100 ng/mL	25	
PCP Phencyclidine	Urine	Strip	DPC-101	✓	✓		25 ng/mL	50	4.0 mm
			DPC-111	✓			25 ng/mL	25 Tests/Canister	
		Device	DPC-102	✓	✓		25 ng/mL	40	
		Dip Card	DPC-114	✓	✓		25 ng/mL	25	
PPX Propoxyphene	Urine	Strip	DPP-101	✓	✓	✓	300 ng/mL	50	4.0 mm
			DPP-111	✓			300 ng/mL	25 Tests/Canister	
		Device	DPP-102	✓	✓	✓	300 ng/mL	40	
		Dip Card	DPP-114			✓	300 ng/mL	25	
TCA Tricyclic Antidepressants	Urine	Strip	DTC-101	✓	✓		1000 ng/mL	50	4.0 mm
			DTC-111	✓			1000 ng/mL	25 Tests/Canister	
		Device	DTC-102	✓	✓		1000 ng/mL	40	
		Dip Card	DTC-114	✓	✓		1000 ng/mL	25	
THC Marijuana	Urine	Strip	DTH-101	✓	✓		50 ng/mL	50	4.0 mm
			DTH-B101	✓			150 ng/mL	50	
			DTH-C101	✓			20 ng/mL	50	
			DTH-111	✓			50 ng/mL	25 Tests/Canister	

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
THC Marijuana	Urine	Device	DTH-102	✓	✓		50 ng/mL	40	4.0 mm
			DTH-B102	✓			150 ng/mL	40	
			DTH-C102	✓			20 ng/mL	40	
		Dip Card	DTH-114	✓	✓		50 ng/mL	25	
			DTH-B114	✓			150 ng/mL	25	
			DTH-C114	✓			20 ng/mL	25	
EDDP Methadone Metabolite	Urine	Strip	DED-101	✓	✓		300 ng/mL	50	4.0 mm
			DED-A101	✓			100 ng/mL		
		Device	DED-102	✓			300 ng/mL	40	
			DED-A102	✓			100 ng/mL		
		Dip Card	DED-114	✓			300 ng/mL	25	
			DED-A114	✓			100 ng/mL	25	
TRA Tramadol	Urine	Strip	DTR-101	✓			100 ng/mL	5	4.0 mm
			DTR-111	✓			100 ng/mL	25	
		Device	DTR-102	✓			100 ng/mL	40	
		Dip Card	DTR-114	✓			100 ng/mL	25	
Multi-Drug 2 Drug Screen Test Panel	Urine	2 Drugs	DOA-124	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 3 Drug Screen Test Panel	Urine	3 Drugs	DOA-134	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 4 Drug Screen Test Panel	Urine	4 Drugs	DOA-144	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 5 Drug Screen Test Panel	Urine	5 Drugs	DOA-154	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 6 Drug Screen Test Panel	Urine	6 Drugs	DOA-164	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 7 Drug Screen Test Panel	Urine	7 Drugs	DOA-174	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 8 Drug Screen Test Panel	Urine	8 Drugs	DOA-184	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 9 Drug Screen Test Panel	Urine	9 Drugs	DOA-194	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 10 Drug Screen Test Panel	Urine	10 Drugs	DOA-1104	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 11 Drug Screen Test Panel	Urine	11 Drugs	DOA-1114	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug 12 Drug Screen Test Panel	Urine	12 Drugs	DOA-1124	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Screen Test Panel and Speciment Validity Test*	Urine	2-12 Drugs	DUA-1X4	✓	✓		See Insert	25	4.0 mm
Multi-Drug Screen Test Panel and Speciment Validity Test*			DUB-1X4	✓	✓		See Insert	25	4.0 mm
Multi-Drug Screen Test Panel and Speciment Validity Test*			DUD-1X4	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 2 Drug Screen Test Device	Urine	2 Drugs	DOA-125	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 3 Drug Screen Test Device	Urine	3 Drugs	DOA-135	✓	✓	✓	See Insert	25	4.0 mm

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
Multi-Drug Multi-Line 4 Drug Screen Test Device	Urine	4 Drugs	DOA-145	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 5 Drug Screen Test Device	Urine	5 Drugs	DOA-155	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 6 Drug Screen Test Device	Urine	6 Drugs	DOA-165	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 7 Drug Screen Test Device	Urine	7 Drugs	DOA-175	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 8 Drug Screen Test Device	Urine	8 Drugs	DOA-185	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 9 Drug Screen Test Device	Urine	9 Drugs	DOA-195	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 10 Drug Screen Test Device	Urine	10 Drugs	DOA-1105	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 11 Drug Screen Test Device	Urine	11 Drugs	DOA-1115	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 12 Drug Screen Test Device	Urine	12 Drugs	DOA-1125	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 2 Drug Screen Test Panel with Integrated Cup	Urine	2 Drugs	DOA-127	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 3 Drug Screen Test Panel with Integrated Cup	Urine	3 Drugs	DOA-137	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 4 Drug Screen Test Panel with Integrated Cup	Urine	4 Drugs	DOA-147	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 5 Drug Screen Test Panel with Integrated Cup	Urine	5 Drugs	DOA-157	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 6 Drug Screen Test Panel with Integrated Cup	Urine	6 Drugs	DOA-167	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 7 Drug Screen Test Panel with Integrated Cup	Urine	7 Drugs	DOA-177	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 8 Drug Screen Test Panel with Integrated Cup	Urine	8 Drugs	DOA-187	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 9 Drug Screen Test Panel with Integrated Cup	Urine	9 Drugs	DOA-197	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 10 Drug Screen Test Panel with Integrated Cup	Urine	10 Drugs	DOA-1107	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 11 Drug Screen Test Panel with Integrated Cup	Urine	11 Drugs	DOA-1117	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 12 Drug Screen Test Panel with Integrated Cup	Urine	12 Drugs	DOA-1127	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with Integrated Cup and Specimen Validity Test*	Urine	2-12 Drugs	DUA-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with Integrated Cup and Specimen Validity Test*			DUB-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with Integrated Cup and Specimen Validity Test*			DUD-1X7	✓	✓	✓	See Insert	25	4.0 mm
Multi-Drug Multi-Line 2 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	2 Drugs	DBO-127	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 3 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	3 Drugs	DBO-137	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 4 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	4 Drugs	DBO-147	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 5 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	5 Drugs	DBO-157	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 6 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	6 Drugs	DBO-167	✓	✓		See Insert	25	4.0 mm

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
Multi-Drug Multi-Line 7 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	7 Drugs	DBO-177	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 8 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	8 Drugs	DBO-187	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 9 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	9 Drugs	DBO-197	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 10 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	10 Drugs	DBO-1107	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 11 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	11 Drugs	DBO-1117	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line 12 Drug Screen Test Panel with E-Z Split Key™ Cup II	Urine	12 Drugs	DBO-1127	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with E-Z Split Key™ Cup II and Specimen Validity Test*	Urine	2-12 Drugs	DUA-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with E-Z Split Key™ Cup II and Specimen Validity Test*			DUB-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Drug Multi-Line Drug Screen test Panel with E-Z Split Key™ Cup II and Specimen Validity Test*			DUD-1X7	✓	✓	✓	See Insert	25	4.0 mm
Multi-Line Screen 2 Test Panel with Integrated iCup®	Urine	2 Drugs	DGO-127	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen 3 Test Panel with Integrated iCup®	Urine	3 Drugs	DGO-137	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen 4 Test Panel with Integrated iCup®	Urine	4 Drugs	DGO-147	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen 5 Test Panel with Integrated iCup®	Urine	5 Drugs	DGO-157	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen 6 Test Panel with Integrated iCup®	Urine	6 Drugs	DGO-167	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 7 Test Panel with Integrated iCup®	Urine	7 Drugs	DGO-177	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 8 Test Panel with Integrated iCup®	Urine	8 Drugs	DGO-187	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 9 Test Panel with Integrated iCup®	Urine	9 Drugs	DGO-197	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 10 Test Panel with Integrated iCup®	Urine	10 Drugs	DGO-1107	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 11 Test Panel with Integrated iCup®	Urine	11 Drugs	DGO-1117	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen 12 Test Panel with Integrated iCup®	Urine	12 Drugs	DGO-1127	✓	✓		See Insert	25	3.5/4.0 mm
Multi-Line Screen Test Panel with Integrated iCup® and Specimen Validity Test*	Urine	2-12 Drugs	DUA-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen Test Panel with Integrated iCup® and Specimen Validity Test*			DUB-1X7	✓	✓		See Insert	25	4.0 mm
Multi-Line Screen Test Panel with Integrated iCup® and Specimen Validity Test*			DUD-1X7	✓	✓	✓	See Insert	25	4.0 mm
E-Z Start Multi-Drug 2 Test Cup	Urine	2 Drugs	DEO-127	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 3 Test Cup	Urine	3 Drugs	DEO-137	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 4 Test Cup	Urine	4 Drugs	DEO-147	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 5 Test Cup	Urine	5 Drugs	DEO-157	✓	✓		See Insert	25	4.0 mm

## DRUGS OF ABUSE (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Cut-off Level	Tests per kit	Dimension
E-Z Start Multi-Drug 6 Test Cup	Urine	6 Drugs	DEO-167	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 7 Test Cup	Urine	7 Drugs	DEO-177	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 8 Test Cup	Urine	8 Drugs	DEO-187	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 9 Test Cup	Urine	9 Drugs	DEO-197	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 10 Test Cup	Urine	10 Drugs	DEO-1107	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 11 Test Cup	Urine	11 Drugs	DEO-1117	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug 12 Test Cup	Urine	12 Drugs	DEO-1127	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug Test Cup and Specimen Validity Test*	Urine	2-12 Drugs	DUA-1X7	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug Test Cup and Specimen Validity Test*			DUB-1X7	✓	✓		See Insert	25	4.0 mm
E-Z Start Multi-Drug Test Cup and Specimen Validity Test*			DUD-1X7	✓	✓	✓	See Insert	25	4.0 mm
Multi-CLIN Drug Screen Test Device	Urine	Device	DCB-1115		✓	✓	See Insert	25	4.0 mm
Adulteration Strips (Creatinine/Nitrite/Glutaraldehyde/pH/Specific Gravity/Oxidant) *	Urine	6 parameters	DUC-111				See Insert	25	4.0 mm

\* Adulteration / Specimen Validity Tests are not regulated by FDA    \*\* Class I specimen transport and storage container (exempt from premarket notification)

## WOMEN'S HEALTH

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
hCG Pregnancy	Urine	Strip	FHC-101	✓	✓	✓	25 mIU/mL	50	4.0 mm
			FHC-U101	✓		✓	10 mIU/mL		
			FHC-111	✓	✓	✓	25 mIU/mL	25*	

## WOMEN'S HEALTH (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
hCG Pregnancy	Urine	Strip	FHC-U111	✓		✓	10 mIU/mL	25*	4.0 mm
			FHC-101H	✓	✓		25 mIU/mL	1/2	
			FHC-U101H	✓			10 mIU/mL		
		Device	FHC-102	✓	✓	✓	25 mIU/mL	40	4.0 mm
			FHC-U102	✓		✓	10 mIU/mL		
			FHC-P102				25 mIU/mL	25	7.2 mm
			FHC-102H	✓	✓		25 mIU/mL	1	4.0 mm
		FHC-U102H	✓			10 mIU/mL			
		Midstream	FHC-103	✓	✓	✓	25 mIU/mL	1/2	7.2 mm
			FHC-U103	✓			10 mIU/mL	1/2	
	FHC-P103		✓	✓		25 mIU/mL	1/2	8.0 mm	
	FHC-T103		✓			20 mIU/mL	1/2	8.0 mm	
	Dip Card	FHC-104				25 mIU/mL	20	4.0 mm	
	FHC-114	✓				25 mIU/mL	25	4.0 mm	
	Urine/Serum	Strip	FHC-201	✓	✓	✓	25 mIU/mL	50	4.0 mm
			FHC-U201	✓			10 mIU/mL	50	
			FHC-211	✓	✓	✓	25 mIU/mL	25*	
			FHC-U211	✓			10 mIU/mL		
Device		FHC-202	✓	✓	✓	25 mIU/mL	40		
		FHC-U202	✓		✓	10 mIU/mL			
LH Ovulation	Urine	Strip	FLH-101	✓	✓	✓	40 mIU/mL	50	4.0 mm
			FLH-U101	✓			30 mIU/mL		
			FLH-U101H	✓			30 mIU/mL	5	
			FLH-111	✓			40 mIU/mL	25*	
		FLH-101H	✓			40 mIU/mL	5		
	Device	FLH-102	✓	✓	✓	40 mIU/mL	40		
		FLH-U102	✓			30 mIU/mL			

## WOMEN'S HEALTH (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
LH Ovulation	Urine	Device	FLH-102H	✓			40 mIU/mL	5	4.0 mm
			FLH-U102H	✓			30 mIU/mL		
		Midstream	FLH-103	✓	✓	✓	40 mIU/mL	5	7.2 mm
			FLH-U103	✓			30 mIU/mL		
FSH Menopause	Urine	Device	FFS-102	✓		✓	25 mIU/mL	40	4.0 mm
			FFS-102H	✓			25 mIU/mL	2	
		Midstream	FFS-103	✓	✓		25 mIU/mL	2	7.2 mm

\*Also available in canisters of 20 or 30

## INFECTIOUS DISEASES

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
Adenovirus	Feces	Device	IAD-602	✓			See insert	25	4.0 mm
Chlamydia	Swab/Urine	Device	ICH-502	✓			See insert	20**	4.0 mm
			ICH-502S	✓			See insert	25	
			ICH-502C				See insert	20	
CMV IgM Device	Serum/Plasma	Device	ICM-302				See insert	40	4.0 mm
Dengue	Whole Blood/ Serum/Plasma	Strip	IDE-401				See insert	20	4.0 mm
		Device	IDE-402	✓			See insert	10	
			IDE-402H				See insert	10	
Gonorrhea	Swab	Strip	IGO-501				See insert	25	4.0 mm
		Device	IGO-502				See insert	20	
Gonorrhea and Chlamydia Combo	Swab	Combo Device	IGC-525				See insert	20	4.0 mm
HBV Hepatitis B Virus Combo	Serum/Plasma	Device	IHB-355				See Insert	25	4.0 mm
HBsAb Hepatitis B Surface Antibody	Serum/Plasma	Strip	IHBsb-301				See insert	50	4.0 mm
			IHBsb-311					25*	
		Device	IHBsb-302					40	
	Whole Blood/ Serum/Plasma	Strip	IHBsb-401					50	3.0 mm

## INFECTIOUS DISEASES (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
HBsAg Hepatitis B Surface Antigen	Serum/Plasma	Strip	IHBsg-301				See insert	50	4.0 mm
			IHBsg-311					25*	
		Device	IHBsg-302					40	
			IHBsg-302C					40	
	Whole Blood/ Serum/Plasma	Strip	IHBsg-401					50	3.3 mm
			IHBsg-411					25*	
	Device	IHBsg-402				40	4.0 mm		
HBsAg Ultra Hepatitis B Surface Antigen	Serum/Plasma	Strip	IHBsg-U301				See insert	50	4.0 mm
			IHBsg-U311					25*	
		Device	IHBsg-U302					40	
HBsAb Hepatitis B Envelope Antibody	Serum/Plasma	Device	IHBeb-302				See insert	40	4.0 mm
HBsAg Hepatitis B Envelope Antigen	Serum/Plasma	Device	IHBeg-302				See insert	40	4.0 mm
HBcAb Hepatitis B Core Antibody	Serum/Plasma	Device	IHBcb-302				See insert	40	4.0 mm
HCV Hepatitis C Virus	Serum/Plasma	Strip	IHC-301				See Insert	50	4.0 mm
			IHC-302					40	
		Device	IHC-302C					40	
	Whole Blood/ Serum/Plasma	Strip	IHC-401					50	
			Device	IHC-402					
HBsAg and HCV Combo	Serum/Plasma	Combo Device	IHBC-325				See insert	40	4.0 mm
HIV 1/2 Human Immunodeficiency Virus 1/3	Serum/Plasma	Device	IHI-U302				See insert	40	4.0 mm
	Whole Blood/ Serum/Plasma	Strip	IHI-401					50	
			IHI-411					25*	
		Device	IHI-402					40	
HIV 1/2/O Tri-line Human Immunodeficiency Virus 1/2/O	Whole Blood/ Serum/Plasma	Device	IHI-T402				See Insert	10 / 40	4.0 mm
			IHI-T402H					10	
			IHI-GT402					10 / 40	
			IHI-GT402H					10	

## INFECTIOUS DISEASES (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
--------------	----------	--------	-------------	-----------	----------------	------------------------	-------------	---------------	-----------

<i>H. pylori</i>	Serum/Plasma	Device	IHP-302	✓		✓	See Insert	40	4.0 mm
	Whole Blood/ Serum/Plasma	Strip	IHP-401	✓	✓	✓		50	
		Device	IHP-402	✓	✓	✓		40	
	Feces			IHP-602	✓				
HSV 1/2 IgM	Serum/Plasma	Device	ISM-302				See Insert	40	4.0 mm
Malaria P.f.	Whole Blood	Strip	IMA-401				See Insert	50	4.0 mm
		Device	IMA-402	✓		✓		40	4.0 mm
Mono	Whole Blood/ Serum/Plasma	Strip	IMO-401		✓		See Insert	20	4.0 mm
			IMO-411					25 Tests/Canister	
		Device	IMO-402	✓	✓			20	
Rubella IgM	Serum/Plasma	Device	IRM-302				See Insert	40	4.0 mm
Rotavirus	Feces	Device	IRO-602	✓			See Insert	25	4.0 mm
Rotavirus and Adenovirus Combo	Feces	Combo Device	IRA-625	✓			See Insert	25	4.0 mm
Strep A	Throat Swab	Strip	IST-501	✓	✓	✓	See Insert	25**	4.0 mm
			IST-A501	✓				10	
			IST-501S	✓				50	
			IST-501H					1**	
			IST-511					25 Tests/Canister	
		Device	IST-502	✓	✓	✓		20**	4.0 mm
			IST-A502	✓				10/25	
			IST-502S	✓		✓		25	
			Twist Device	IST-502T	✓	✓		✓	
Syphilis	Serum/Plasma	Device	ISY-U302	✓			See Insert	40	4.0 mm
	Whole Blood/ Serum/Plasma	Strip	ISY-U401	✓				50	
			ISY-U411					25*	
	Device	ISY-U402	✓			40			

## INFECTIOUS DISEASES (Continued)

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
TorCH IgM	Serum/Plasma	Combo Device	ITOM-345				See Insert	25	4.0 mm
Toxo IgM	Serum/Plasma	Device	ITM-302				See Insert	40	4.0 mm
Tuberculosis	Whole Blood/ Serum/Plasma	Strip	ITB-401	✓			See Insert	50	3.5 mm
		Device	ITB-402	✓				40	4.0 mm

\*Also available in canisters of 20 or 30

\*\*Includes swabs and/or controls

## ONCOLOGY

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
AFP Alpha-Fetoprotein	Whole Blood/ Serum/Plasma	Strip	TAF-401				10 ng/mL	50	4.0 mm
		Device	TAF-402	✓			10 ng/mL	40	4.0 mm
CEA Carcinoembryonic Antigen	Whole Blood/ Serum/Plasma	Strip	TCE-401				5 ng/mL	50	4.0 mm
		Device	TCE-402				5 ng/mL	40	4.0 mm
FOB Fecal Occult Blood	Feces	Strip	TFO-601	✓		✓	50 ng/mL 6 µg/g	25	4.0 mm
			TFO-601S				50 ng/mL 6 µg/g	40	4.0 mm
		Device	TFO-602	✓			✓	50 ng/mL 6 µg/g	25
PSA Prostate Specific Antigen Semi-Quantitative	Whole Blood/ Serum/Plasma	Strip	TPS-401				4 ng/mL (R=10 ng/mL)	50	4.0 mm
		Device	TPS-402					40	4.0 mm
			TPS-U402				3 ng/mL (R=10 ng/mL)	40	4.0 mm
PSA Prostate Specific Antigen Qualitative	Whole Blood/ Serum/Plasma	Device	TPS-UB402				3 ng/mL	40	4.0 mm

## CARDIOLOGY

Product Name	Specimen	Format	Catalog No.	CE Marked	US FDA Cleared	Health Canada Licensed	Sensitivity	Tests per kit	Dimension
CK-MB	Whole Blood/ Serum/Plasma	Device	CCK-402	✓			5 ng/mL	20	4.0 mm
MYO Myoglobin	Whole Blood/ Serum/Plasma	Device	CMY-402	✓			50 ng/mL	20	4.0 mm
cTnl Troponin I	Whole Blood/ Serum/Plasma	Device	CTI-402				0.5 ng/mL	20	4.0 mm
Myoglobin/CK-MB/Troponin I Combo	Whole Blood/ Serum/Plasma	Device	CMA-435	✓			See Insert	20	4.0 mm
CRP C-Reactive Protein	Whole Blood/ Serum/Plasma	Device	CCR-402				1 mg/L	20	4.0 mm
NT-proBNP Rapid Test Device	Whole Blood/ Serum/Plasma	Device	CBN-402				450 pg/L	20	4.0 mm

Telephone: +1 858.805.8900 Fax: +1 858.695.8796

9975 Summers Ridge Road. San Diego, CA 92121 USA

info@innovaconinc.com